Invia® Liberty™

NEGATIVE PRESSURE WOUND THERAPY SYSTEM

- EN** Clinician instructions for use
- DE** Gebrauchsanweisung für medizinisches Personal
- FR** Mode d'emploi à l'usage du personnel médical
- ES** Instrucciones de uso para profesionales sanitarios
- PT** Instruções de utilização para pessoal médico
- NL** Gebruiksaanwijzing voor artsen
- SV** Bruksanvisning för vårdpersonal
- IT** Istruzioni per l'uso per il medico
- FI** Käyttöohjeet lääkärille
- DA** Brugsanvisning til sundhedspersonalet
- NO** Bruksanvisning for helsepersonell

NPWT University

Medela.com website provides all the educational tools you need to set-up and operate the Invia Liberty NPWT system. The Medela University page includes instructions for use, quick cards and training videos to support you in becoming a confident and knowledgeable NPWT user.

www.medela.com

NPWT University can be found at www.medela.com via the “University” tab. It is available on all formats including desktop, smartphone and tablet. We recommend saving NPWT University as a shortcut on desktop, smartphone and tablets for an even quicker access.

Table of contents

Introduction	6
To consider before use	6
Intended use	6
Indications for use	7
Contraindications	7
Warnings, cautions and safety instructions	7
WARNINGS	8
CAUTIONS	9
Safety instructions	10
Wound assessment	10
Safety-related checks	11
Dressing technique	11
Invia Liberty pump	12
Display	13
Invia Liberty pump disposables	14
Invia Liberty canister with solidifier 300ml and 800ml	14
Invia Liberty canister tubing/Y-connector/Drain adapter	15
Power supply	16
USB port	16
Preparation for use	17
Invia Liberty pump user modes	18
Administrative mode	18
Patient mode	18
Administrative Mode	19
Device setup	19
Change pressure level	20
Change pressure set-up mode	21
Select air leakage volume	21
Change therapy mode	23
Change settings	24
Therapy Log file	25
To open the Therapy Log file	25

Patient mode	25
Turn ON	25
Check pressure	25
Air leakage indicator	25
Standby	26
Turn OFF	26
Set up carrying case	26
Change Invia Liberty canister and Invia Liberty tubing	27
Battery charging	29
Alarms	30
“Warning”	31
“Alarm”	31
“Internal fault”	31
Alarm table	32
Accessories overview	36
Wound dressings	37
Sterility and requirements for usage	37
Cleaning and disinfection	38
Disposal	39
Maintenance/Safety-related check	39
Guarantee	39
Service life	39
Signs and symbols	40
Technical specifications	41
Electromagnetic compatibility (EMC)	42
Electromagnetic emissions	42
Electromagnetic immunity	43
Electromagnetic immunity	44
Recommended separation distance	45

Introduction

With the Invia Liberty you have selected a system for use in Negative Pressure Wound Therapy (NPWT). The lightweight Invia Liberty pump provides an adjustable negative pressure range and two therapy modes along with an electronic measuring and monitoring system. The pump is quiet during operation and has optical and acoustic status alarms for patient safety.

Invia Liberty is portable and can be operated independent of the electrical outlet thanks to a rechargeable battery. Acoustic and optical signals are triggered for variances from the set values as well as for faults.

To consider before use

These instructions for use are a general guide for the use of the Invia Liberty pump with associated products. Medical matters must be addressed by a physician.

In order to ensure safe and proper operation of Medela products, a quality management system is used. Please meet the conditions below (failing to do so will void the warranty). Invia Liberty NPWT system is to be used exclusively as described in these instructions for use.

- Before initiating NPWT treatment, read the instructions for use, indications, contraindications, warnings, cautions and safety instructions. Nonobservance and incorrect use can lead to considerable dangers and cause pain and injury to the patient.
- Safe and effective operation of this device requires specific instructions from a physician.
- For use only by persons who have been adequately trained in wound care and negative pressure wound therapy.
- Therapy changes (pressure level, constant or intermittent mode) should only be done as prescribed by a physician.
- In these instructions for use "pressure" in general implies "negative pressure".
- Please keep in mind that each wound is unique and must be assessed by a qualified medical professional who must use his / her best clinical judgment when applying this therapy. The pressure level and therapy mode must be adapted to each individual patient according to his / her medical knowledge and according to the wound healing phase.

Intended use

Intended user

The Invia Liberty NPWT system is intended to be used by healthcare professionals or adequately trained lay users.

Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.

Intended patient population

The Invia Liberty NPWT system is intended to be used on patients only exhibiting conditions as described in the indications for use. The device has not been studied in pediatric patients.

Intended environment

The Invia Liberty NPWT system is intended for use in acute, extended and home care settings.

Indications for use

The Invia Liberty Negative Pressure Wound Therapy (NPWT) system is indicated for patients who would benefit from a suction device (NPWT) as it creates an environment that promotes wound healing by secondary or tertiary (delayed primary) intention by preparing the wound bed for closure, reducing edema, promoting granulation tissue formation and perfusion, and by removing exudate and infectious material. When used on closed surgical incisions, the Invia Liberty NPWT system is also intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of Negative Pressure Wound Therapy.

The Invia Liberty NPWT system is appropriate for the following indications:

- Acute or subacute wounds
- Chronic wounds
- Dehisced wounds
- Pressure ulcers
- Diabetic/Neuropathic ulcers
- Venous insufficiency ulcers
- Traumatic wounds
- Partial thickness burns
- Flaps and grafts
- Closed surgical incisions

Contraindications

The Invia Liberty NPWT system is contraindicated in the presence of:

- Necrotic tissue with eschar present
- Untreated osteomyelitis
- Non-enteric and unexplored fistulas
- Malignancy in the wound
- Exposed vasculature
- Exposed nerves
- Exposed anastomotic site of blood vessels or bypasses
- Exposed organs

Warnings, cautions and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

The Invia Liberty Negative Pressure Wound Therapy system (Invia Liberty NPWT system) is intended for the use described in these instructions for use.

Medela is only responsible for the effect on BASIC SAFETY, reliability and performance of the Invia Liberty NPWT system if it is used in accordance with the instructions for use.

Please read and observe these warnings and safety instructions before operation.

These instructions for use must be kept with the device.

Please note that these instructions for use are a general guide for the use of the product. Medical situations must be addressed by a physician.

WARNINGS

- Do not modify this equipment without authorization from the manufacturer.
- This manual provides general guidelines for the use of the Invia NPWT system.
- The safe and effective operation of this device requires specific instruction from a physician.
- No modification of this equipment is allowed.
- For use only by healthcare professionals who have been adequately trained in suction procedures, wound care, negative pressure wound therapy and in the use of aspirators or adequately trained lay users. Healthcare professionals are responsible to train lay users according to the patient instructions for use and explain all related safety information.
Caution: Incorrect use can cause pain and injury to the patient.
- Consult the indications for use, cautions and contraindications when using the Invia Liberty as a vacuum source with the Invia NPWT system. Failure to read and follow all instructions in this manual prior to use may result in death or injury of the patient.
- Failure to obtain consent and any additional instructions from the treating physician prior to use, may lead to death or injury of the patient.
- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.
- The device is not for use while bathing or showering
- Do not dry Invia Liberty in a microwave.
- Data transfer via USB is not possible in the running mode.
- A patient undergoing NPWT requires frequent supervision. Objective indications or signs of possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent drainage). Monitor the device, wound, surrounding skin and patient status and comfort level frequently to ensure efficient, safe treatment and patient comfort.
- Do not place the foam/gauze dressing directly on exposed blood vessels, organs, nerves, tendons, bones or ligaments. When using the Invia Liberty NPWT system in close proximity to these structures a protective barrier, such as a non-adherent wound contact layer, must be used.
- Serious or fatal injury can result from bone fragments or sharp edges (e.g. staples or hardware) that could puncture protective barriers, vessels or organs.
- Patient must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately stop use of the pump, apply pressure on wound dressing and seek immediate Emergency Medical Attention.
- Should a spinal cord injury patient experience autonomic hyperreflexia, discontinue treatment with the Invia Liberty NPWT system and consult a physician immediately.
- Never place the Invia Liberty pump in water or liquids. Clamp the drain and disconnect from the dressing prior to bathing or showering.
- Consider the use of a protective barrier on skin that may come in contact with the tubing, especially in patients with fragile skin.
- Invia NPWT instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.
- This device has not been studied in pediatric patients.
- Clamp the drain and disconnect the Invia Liberty pump prior to patient entering hyperbaric oxygen chamber (HBO) or Positron Emission Tomography (PET).
- The Invia Liberty NPWT system is not for use in the Magnetic Resonance (MR) environment, so do not take the Invia Liberty NPWT system into this environment.

- Explosion hazard – the Invia Liberty pump is not for use in potentially explosive environments including oxygen enriched environments and in areas of flammable anesthetics.
- In the event that defibrillation is required, disconnect the pump from the wound dressing before the patient is defibrillated.
- Special care is advised for dressing placement and removal in order to avoid situations such as unintentional gauze or foam retention.
- This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100°C, incinerate or dispose of in fire.

Contact your local Medela customer service representative for assistance with product operations.

CAUTIONS

The following statements describe medical conditions that may require special care for the safe and effective use of the Invia NPWT system.

- Patients at high risk for bleeding and hemorrhage.
- Patients taking anticoagulants or platelet aggregation inhibitors, or with patients experiencing active bleeding or difficult wound hemostasis.
- Wounds that involve a fistula.
- Using Invia NPWT system in close proximity to blood vessels and organs or exposed organs, vessels, nerves, tendons, or ligaments. Provide necessary protection of all vessels and organs using a protective barrier.
- Patients with a history of vascular anastomosis or friable, irradiated, sutured or infected blood vessels.
- Use near vagus nerve (bradycardia) or use on patients with a history of spinal cord injury (stimulation of sympathetic nervous system).
- Circumferential dressing application.

Incorrect use may cause pain and injury to the patient. Excessive negative pressure, a too tight adhesive cover dressing, or an infection of the wound may cause pain to the patient. In either case, the dressing must be removed and the wound assessed.

The patient should be monitored regularly according to physician instructions and facility guidelines to monitor patient comfort, therapy compliance and signs of wound infection.

Do not use an Invia Liberty canister or tubing if the sterile packaging is damaged.

The Negative Pressure Wound Therapy must be used 24 hours per day without interruption. If the pump is stopped for more than two hours, the dressing must be changed and therapy restarted.

Consider the patient's size and weight when prescribing this device.

Consider mode of therapy – intermittent versus continuous.

Safety instructions

- Before you plug in the device, please check that your local power supply is the same as the voltage given on the specification plate.
- The Invia Liberty pump is verified within the scope of conformity evaluation and is only to be used with products included in the Invia Liberty NPWT system and distributed by Medela. Medela can only guarantee the effective performance of the system with these products.
- Wireless communication equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, and walkie-talkies can interfere with the Invia Liberty pump and should be kept at a minimum distance of 30 cm away from the Invia Liberty pump.
- The patient should be regularly monitored according to facility or institution guidelines.
- Invia Liberty pump must remain in an upright position during use.
- Supervision is necessary when the Liberty pump is used in the vicinity of children.
- Do not use Invia Liberty pump if:
 - The power cord or plug are damaged
 - The device is not functioning properly
 - The device is damaged
 - The device has apparent safety defects
- Never pull the plug out of the main socket by pulling on the connecting cable.
- Keep the Invia Liberty pump with associated products away from hot surfaces.
- Never place the Invia Liberty pump, charger or docking station device in water or other liquids and keep the charger connector away from moisture or immersion in water.
- The Invia Liberty pump must not be used for suctioning explosive, easily flammable or corrosive liquids.
- The tubing connected to the canister must never come in direct contact with the wound area.

Wound assessment

CAUTION

Patient monitoring: The patient should be monitored regularly according to the physician's instructions and facility guidelines to check for patient comfort, therapy compliance and signs of infection.

WARNING

Objective indications or signs of a possible infection or complication must be addressed immediately (e.g. fever, pain, redness, increased warmth, swelling or purulent discharge). Non-observation can lead to considerable danger to the patient.

Observe wound/periwound tissue and exudate for signs of infection or other complications. Most common signs of infection include redness, tenderness, fever, swelling, itching, increased warmth in the wound area, strong odor or purulent discharge. Additional symptoms include nausea, vomiting, diarrhea, headache, dizziness, fainting, sore throat with swelling of the mucous membranes, disorientation, high fever (>38.8° C), refractory hypotension, orthostatic hypotension, or erythroedema (a sunburn-like rash). More serious complications of infection include pain, discomfort, fever, gangrene, toxic or septic shock. If more serious complications of infection occur, discontinue therapy and consult a physician immediately.

Safety-related checks

For the safety-related checks, the device should be maintained and repaired throughout its service life in compliance with the service procedures.

The Invia Liberty pump is a device in protection class II (EN IEC 60601-1), the safety-related checks are confined to visual inspection of the housing and charger for damage. This check must be performed prior to each use.

Devices of protection class II do not have a protective earth conductor; there is therefore no need to check the earth leakage current.

The Invia Liberty pump enclosure is made entirely of insulated material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

Even when suctioning a conductive fluid until the overflow protection device activates, measurements of the patient leakage current using common measuring instruments will not reveal measurable values.

Invia Liberty pump does not have patient circuits or functional earth connections.

Dressing technique

Consult the appropriate Invia wound dressing instructions for use for information regarding dressing applications. Perform a thorough wound cleansing per physician orders or facility protocol prior to dressing applications.

Dressing change



 For information regarding the dressing change, consult the specific Invia wound Dressing Instructions for Use.

When treating infected wounds or wounds more susceptible to tissue in-growth into wound filler material, more frequent dressing changes may be needed. The frequency of dressing changes should be based on an evaluation of the wound characteristics rather than standard recommendations.


 Important things to remember:

- Routinely check that the negative pressure level and the therapy mode is at the prescribed setting.
- Monitor the integrity of the dressing and that the pump is on and running.

Invia Liberty pump


The max flow of Invia Liberty pump is 5 liters/minute with an adjustable pressure range of -40 to -200 mmHg (-5 to -27 kPa).

Display

Invia Liberty pump disposables

Invia Liberty canister with solidifier 300 ml and 800 ml

Material:

Polypropylene

Accuracy of graduation:

+/- 2.5 % (in the upright position)

Composition of solidifier:

Cross-linked sodium polyacrylate

CAUTIONS


For appropriate, safe operation the Invia Liberty pump must remain in an upright position during use.

If the pump tips over, set it upright again. The special construction of the safety chamber in the upper region of the canister protects the overflow protection/bacteria filter from clogging immediately, if tipped over.

When the canister is full and the pump tips over, this function is made inoperative, since the secretions will flow into the safety chamber and clog up the hydrophilic filter. In this case an alarm will sound and the canister will need to be replaced.

Invia Liberty canister tubing/Y-connector/Drain adapter

EN

Material: TPE
Length: 1.5 m
Diameter: 2.0/3.0 mm

The canister tubing is comprised of two lumens:

The smaller lumen (measuring tubing) regulates the pressure and the larger lumen (suction tubing) removes the fluid from the wound into the canister. A hydrophilic overflow protection/bacteria filter in the tubing base helps to prevent contamination of Invia Liberty pump.

The Invia Y-connector with Quick-connector

For warnings and cautions of the Invia Y-connector with Quick-connector consult the Invia Y-connector with Quick-connector instructions for use.

Power supply

CAUTIONS

Before you charge the device, please check that your local power supply is the same as the voltage given on the specification plate.

WARNINGS

This product contains lithium-ion batteries which bear risk of fire, explosion and burns.
Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.

Invia Liberty pump is operable while connected to the electrical power supply or by an internal rechargeable lithium-ion battery. While in use and connected to the electrical power supply, the battery is recharged.

The charge on the battery is dependent upon the run-time of the pump. The run time refers to the effective operation of the motor. Invia Liberty pump motor only turns on when the measured pressure is lower than the set pressure.

Under typical operating conditions, the battery run-time generally exceeds 14 hours. Battery run-time is influenced by the size of the wound, the air leakage in the system and the set pressure. If there is a leakage in the system, the pump motor will run more often, which will reduce the battery run-time.

USB port

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (e.g. IEC 60950 for data processing equipment).

Furthermore all configurations shall comply with the requirements for medical electrical systems (see IEC 60601-1-1 or clause 16 of the 3Ed. of IEC 60601-1, respectively). Anyone connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local representative or the Medela customer service department.


Preparation for use

WARNING

Use only after instruction by trained personnel. Wear gloves for all operations and utilize universal precautions.

1. Check necessary parts

2. Connect Invia Liberty canister tubing

A. Open the tubing packaging


B. Insert the tubing base into the pump as shown (straight push)


3. Click in Invia Liberty canister


A. Unpack new canister and remove the safety sticker. Slip the pegs at the bottom of the canister into the slots at the base of the pump. Lift the canister towards the pump until it clicks into place.

4. Connect the dressing tubing to the canister tubing.

5. Switch Invia Liberty pump ON by pressing []

A. The pump will do a self-test and start running immediately. You should hear the motor running for a short time until the pressure is built up. If the motor keeps running for more than 30 sec. – check system and retry.

 If the above check was successful, switch OFF Invia Liberty pump by pressing [ > 3 seconds]. If the self-test is not successful, follow instructions shown on the display or see the Alarms chapter.

For dressing applications, please refer to the instructions for use provided with the Invia dressings. Turn the Invia Liberty pump on as described in the user modes section.

Invia Liberty pump user modes

Administrative mode

Used by healthcare professional to either set up new patient therapy or to change pump settings such as pressure, air leakage volume, therapy Constant and Intermittent modes, intermittent times, language and time zone. You can enter the Administrative mode when first turning the pump on or while therapy is running.

Patient mode

In Patient mode the pump can be turned on and off, placed in standby mode and an alarm can be muted. In this mode the canister and tubing can be changed and the pump can be charged. When the pump is turned off in this mode, last settings are used by default.

Administrative Mode

 The factory default settings of the Invia Liberty pump are -125 mmHg and constant mode.


CAUTION

The pressure level should always be set according to prescribing healthcare professional instruction.

When an Invia FitPad is used, the set pressure is controlled at the wound site.

Device setup

New patient therapy

1. Switch Invia Liberty pump ON in Administrative mode, press and hold [], and then press [].

Self-test starts

If the self-test is not successful, follow instructions shown on the display or see the Alarm chapter.

2. Acknowledge disclaimer. Press “OK” [] to confirm.

Disclaimer	
This device may only be used by a medically trained person.	
Please confirm that you have read the instructions for use.	
OK	

“New patient? Yes/No”

“Yes” means that Invia Liberty pump will issue a new therapy ID number (= new patient).

This number can be noted in the patient’s file.

The therapy ID number is displayed in the information field.

If yes

Press “Yes” [] to confirm.

Press “OK” [] to enter the main display.


New patient?	
No	Yes
OK	

“No” means that the therapy ID number and settings remains unchanged (= same patient).

Please write the numbers below on Patient file:	
Serial number:	152887
Therapy ID:	1
OK	


If no

Press “No” [] to confirm and enter the main display.


3. Press the Selection buttons [] or [] to set pressure level.
4. Press "OK" [] to confirm and enter the main display.

Administrative mode

1. To enter the Administrative mode, unlock the display, press and hold [], and then press [].

Change pressure level


CAUTIONS

Set the pressure level and therapy mode according to prescribing healthcare professional instruction.

Default pressure mode is Standard with preset pressure levels:

mmHg	-40	-60	-80	-100	-125	-150	-175	-200
(kPa)	-5	-8	-11	-13	-17	-20	-23	-27

1. To enter the Administrative mode, unlock the display, press and hold [], and then press [].
2. Press the Selection buttons [] or [] to set pressure level.
3. Press "OK" [] to confirm.
4. To return to the main display, press "Back" [].

CAUTIONS

If the pressure level is not confirmed, the pump will switch back to previous settings and return automatically to the main display after 5 seconds.

Change pressure set-up mode

Detailed pressure mode allows to select pressure level from -40 to -200 mmHg in increments of 5mmHg (5kPa to 27kPa in increments of 1kPa)


Settings 1/3	
Pressure	-125
Unit Pressure	mmHg
Pressure set-up	Standard
Air Leakage	Standard
Language	English
Back	OK

Settings 1/3	
Pressure	-125
Unit Pressure	mmHg
Pressure set-up	Detailed
Air Leakage	Standard
Language	English
Back	OK

- To enter the Administrative mode, unlock the display, press and hold [], and then press [].
- Press the soft buttons [ + ] at the same time to enter the Settings menu.
- To switch pressure modes, scroll down with [] and select Pressure set-up.
- Press "OK" [] to select.
- To change between Standard and Detailed pressure modes, press Selection buttons [] or [] to choose Standard or Detailed mode.
- Press "OK" [] to confirm.
- Refer to Change pressure level instructions to select pressure values.

Select air leakage volume

 High air leakage mode is recommended when X-large dressing is used.


Default air leakage volume is Standard. If the vacuum on the wound is < 75% of the set pressure level or the air flow is > 1500ml/min, leakage indicator is "full" and flashes, indicating there is a big air leak in the system. The air leakage alarm will go off within 2 minutes if the air flow is above the defined level or after 1 minute if the set pressure is not maintained. Follow the instructions shown on the display, or see the Alarms chapter.


For wounds with air leakage volume between 1500ml/min and 2200ml/min, High air leakage volume mode is available. If the vacuum on the wound is < 75% of the set pressure level or the air flow is > 2200ml/min value, leakage indicator is "full" and flashes, indicating there is a big air leak in the system. The air leakage alarm will go off within 2 minutes if the air flow is above the defined level (2200ml/min) or after 1 minute if the set pressure is not maintained. There will be no leakage alarm during the first 5 minutes which allows establishing the set pressure on a wound with high volume without an alarm.

Selecting between Standard and High air leakage volume:


1. To enter the Administrative mode, unlock the display, press and hold [], and press [].
2. Press the soft buttons [+] at the same time to enter the Settings menu.
3. To change air leakage volume, scroll down with [] and select Air Leakage.
4. Press “OK” [] to select.


5. To change between Standard or High air leakage, choose the mode with the Selection buttons [] or [].

6. Press “OK” [] to confirm.
7. To return to the main display, press “Back” [].


Change therapy mode

- To enter the Administrative mode, unlock the display, press and hold [▶], and then press [].
- To change the therapy mode, press "Standby" [◀ 3 seconds] to put the Invia Liberty pump into the Standby mode.

Constant mode – C


The default therapy mode is Constant.


- To change from Constant mode to Intermittent mode, press [▶] "Change to Intermittent" and then press "On" [◀].

Intermittent mode – I

The default times for Intermittent mode are 5 minutes on and 2 minutes off.

- To change from Intermittent mode to Constant mode, press [▶] "Change to Constant" and press "On" [◀].


1 minute after the last button has been pressed, the Invia Liberty pump switches into Patient mode automatically and the display is locked.

1 minute after the last button has been pressed, the backlight turns off. When pump alarms or any button is pressed the backlight will illuminate.

The settings can only be changed when the pump is in the Administrative mode.

The settings possible to change are, Unit pressure, Air leakage volume, Constant and Intermittent modes, intermittent times, language and time zone. Pump number, pump run-time and pump version can only be viewed but not changed.

Change settings

- To enter the Administrative mode, unlock the display, press and hold [], and then press [].
- Press the soft buttons [+] at the same time to enter the Settings menu.
- To change a setting, choose with the Selection buttons [] or [].

Settings 1/3	
Pressure	-125
Unit Pressure	mmHg
Pressure set-up	Standard
Air Leakage	Standard
Language	English
Back	OK

Selection
-40 to -200 mmHg (-5 to -27 kPa)
kPa/mmHg
Standard / Detailed
Standard / High
Languages according to list

Settings 2/3	
On time	1 min
Off time	1 min
Time zone	+1 h
Pump Initiative	35-28%
Pump Initiative	27 h
Back	OK

Selection
1 – 8 min
1 – 8 min
+12 until -12 hours GMT
as displayed
as displayed

Settings 3/3	
Warnnn	177
Back	OK



Selection
as displayed

- Press “OK” [] to select.
- To change value, choose with the Selection buttons [] or [].
- Press “OK” [] to confirm.

If the changed setting is not confirmed with “OK”, the pump will switch back to the previous setting.

- Press “Back” [] to exit settings.


1 minute after the last button has been pressed, the Invia Liberty pump switches into Patient mode automatically and the display is locked.


Therapy Log file

In the Therapy Log file information regarding runtimes (on/off), pressure settings, therapy modes, alarms, warnings and errors is listed. The last 51 events are logged and displayed.

To open the Therapy Log file

List


Graph

1. To enter the Administrative mode, unlock the display, press and hold [], and then press [].
2. Press Selection buttons [] + [] simultaneously.
3. To view additional pages, scroll with the Selection buttons [] or [].
4. To view the Therapy Log file as a graph, press "Graph" [].
5. To exit the Therapy Log file, press "Back" [].

The Invia Liberty pump switches automatically into main display (in Administrative mode) 30 seconds after the last button has been pressed. After additional 30 seconds, the pump switches into Patient mode and the display is locked.

Patient mode

Turn ON

Switch Invia Liberty pump ON by pressing []


Self-test starts

If the self-test is not successful, follow instructions shown on the display or see the Alarms chapter.

When the pump is turned on in this mode, last settings (therapy mode and pressure level) are used by default.

Check pressure


The set pressure will be shown on the display. The motor will run for a few seconds to build up the pressure. If it runs continuously for more than 30 seconds, check system for leaks.


Air leakage indicator

An air leakage indicator is shown on the display to visualize if there is an air leakage in the system.

a) When the indicator is "empty", the system is air tight.

- b) When the indicator is "half full" there is an air leak in the system, but the pressure and therapy is maintained in accordance with the set pressure.
- c) When the indicator is "full" and flashes, there is a big air leak in the system. When pump is first turned on and a big air leak is present the air leakage alarm will go off within 2 minutes in the Standard air leakage mode or within 5 minutes in the High air leakage mode if the set pressure is not maintained. When therapy is running and a big air leak is detected, alarm will go off within 2 minutes. Follow the instructions shown on the display, or see the Alarms chapter.

Standby

Change the Invia Liberty pump into the Standby mode.

Press "Standby" [> 3 seconds] and the pump motor will stop running.

If the pump is in Standby mode for more than 5 minutes, a warning will go off, follow instructions shown on the display or see the Alarms chapter.

Turn OFF

Press [> 3 seconds] and the pump will turn off.

Set up carrying case

The Invia Liberty is intended to be worn in a carrying case during operation.

1. Open the carrying case.
2. Slide the Invia Liberty into the designated pocket.
3. Close the velcro over the pump.

The shoulder strap can also be used as a belt strap.

To accomplish this, release the strap at the sides of the carrying case and pull it through the two loops located on the back of the carrying case.

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Change Invia Liberty canister minimum once a week, when canister looks full on visual inspection or when "Canister full" / "Filter clogged" alarm signal is activated.


Change Invia Liberty tubing minimum once a week, when tubing looks clogged on visual inspection or when "System clogged" alarm signal is activated.

1. Provide sterile canister and sterile tubing.
2. Press "Standby" [ > 3 seconds] for the pump motor to stop running. Clamp the canister tubing.
 - 3.1 Position the clamp next to the quick-connector and close the clamp on the pump tubing and close the clamp on the dressing tubing.

- 3.2 Disconnect the canister tubing from dressing tubing.


- 3.3 Protect the end of dressing tubing from contamination.

4. Release and remove canister.


5. Seal used canister with a cap located on the side of the canister.

6. Remove canister tubing in direction of the arrow.


7. Unpack new canister tubing and connect to Invia Liberty pump. Insert the tubing base into the pump (straight push).

- 8.1 Unpack new canister and remove the safety sticker.
 - 8.2 Slip the pegs at the bottom of the canister into the slots at the base of the pump. Lift the canister towards pump until it clicks into place.
 - 8.3 Connect the canister tubing to dressing tubing. **Ensure that the tubing is un-clamped.**
 - 8.4 Press "ON" with []. Pressure will begin to built up.
9. Dispose canister and canister tubing in accordance with local procedures. In the home care setting, return disposables to care giver for correct disposal.

Battery charging

The battery must be fully charged before first usage. To charge the battery, plug the charger into the electrical outlet port in the Invia Liberty pump. You can continue to use the pump while it is charging.

The capacity of the battery is not negatively affected if it is charged when partially discharged and the pump does not need to be fully charged at each occasion.

The battery can be charged when appropriate, or in accordance with alarm signal for "Battery low" or "Battery fully discharged (empty)" and instructions in the pump display.

There is approximately 30 minutes of charge remaining on the battery at the onset of the "Battery low" warning. If the battery is completely discharged it will take 3 to 4 hours to charge to full status.

If the Invia Liberty pump is fully charged AND the pump is still connected to an electrical source, will appear in the battery icon.

If the pump is disconnected from an electrical source, 4 bars are visible in the battery icon which indicates that the battery is fully charged.

If Invia Liberty pump has not been in use, the battery must be charged approximately once every 6 months to ensure optimum function.

To disconnect Invia Liberty pump from the electrical source, remove the charger connector and close the safety cover.


WARNING

Do not pull on the cable or the anti-bend protection.

Alarms

The pump distinguishes between “Warning”, “Alarm” and “Internal fault”. If the Invia Liberty pump detects any situations where the therapy cannot be maintained, an acoustic alarm sounds, a fault number and a description of the problem appears on the display. For explanation of the fault number, see the Alarm Table in this chapter.

Example:


CAUTIONS


Invia Liberty instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.

“Warning”

An acoustic alarm sounds and the fault number is shown on the display.

**“Alarm”**

An acoustic alarm sounds and the fault number is shown on the display.

When a Warning/Alarm goes off, an acoustic alarm sounds. A description of the “Warning” or “Alarm” will be shown on the display.

1. Press “Mute” [] to mute and acknowledge the alarm.
The acoustic alarm will resume in 1 minute if the problem is not solved.
2. Follow the instruction shown on the display or see the Alarm table.
3. If the problem cannot be solved, turn off the Invia Liberty pump [ > 3 seconds] and contact the Medela customer service for further instructions.















“Internal fault”

Pump operation stops and an acoustic alarm sounds, “Internal fault” is shown on the display.











1. Press [ > 3 seconds] and the pump will be turned off.
2. Restart the pump by pressing [] and the pump will be turned on.
3. If the Internal fault alarm remains, turn off [ > 3 seconds] Invia Liberty pump and contact the Medela customer service.

Alarm table

	Fault number	Problem description on the display	Troubleshooting on the display	Pressure
Warning	401	 <p>Battery low</p> <p>Mute Standby (3 sec.)</p> <p>401</p>	 <p>Charge battery</p> <p>Mute Standby (3 sec.)</p>	✓
	402	 <p>USB connection not permitted</p> <p>Mute Standby (3 sec.)</p> <p>402</p>	 <p>Unplug USB cable</p> <p>Mute Standby (3 sec.)</p>	
	405	 <p>Standby mode</p> <p>Mute On</p> <p>405</p>	 <p>Switch pump on or off</p> <p>On</p>	✗
	406	 <p>Internal temperature high</p> <p>Mute Standby (3 sec.)</p> <p>406</p>	 <p>Remove the pump from the heat source (e.g. direct sunlight) or remove any additional coverage (e.g. blanket).</p> <p>Mute Standby (3 sec.)</p>	✓
Alarm	301	 <p>Air leak in system</p> <p>Mute Standby (3 sec.)</p> <p>301</p>	 <p>Check dressing for air leakage and if canister is properly inserted.</p> <p>Consult IHU for further instructions.</p> <p>Mute Standby (3 sec.)</p>	✓
	302	 <p>System clogged</p> <p>Mute Standby (3 sec.)</p> <p>302</p>	 <p>Check that tubing is clear and not kinked. Check if canister is full.</p> <p>Consult IHU for further instructions.</p> <p>Mute Standby (3 sec.)</p>	✓
	305	 <p>Battery empty</p> <p>Mute</p> <p>305</p>	 <p>Charge battery</p> <p>Mute</p>	✗
	306	 <p>Canister full</p> <p>Mute Standby (3 sec.)</p> <p>306</p>	 <p>Change canister</p> <p>Mute Standby (3 sec.)</p>	✓

	Remarks/potential cause of fault
	<p>Recharge the battery either by placing the Invia Liberty pump in the docking station or plug in the charger to the electrical outlet port on the pump.</p> <p>Remaining time of battery is approximately 30 minutes.</p>
	<p>Unplug USB cable</p>
	<p>If the pump is in Standby mode for more than 5 minutes, an alarm will go off.</p> <p>To continue therapy press “On” [] or switch off the pump by pressing [ > 3 seconds].</p>
	<p>Cool Invia Liberty pump down as per instructions on the display.</p>
	<p>Dressing: - Check dressing for air leakage. Press firmly around the edges of the dressing, around the drain tubing or on the Invia FitPad. - Apply some additional film dressing to seal the leaking area.</p> <p>Connectors: - Ensure that the dressing tubing is connected properly to the canister tubing. - Ensure that the canister tubing is inserted straight into the pump.</p> <p>Canister: - Ensure that the canister is properly inserted, release the canister and reposition. - Ensure that the orange O-ring/gasket, placed beside the canister tubing on the pump is not missing. Additional O-ring is available via Medela customer service.</p>
	<p>Tubing: - Ensure that the tubing is not twisted, kinked or clamped. - If the canister tubing is clogged, change the tubing.</p> <p>Canister: - If canister is full or filter clogged, replace canister. - If dressing tubing is clogged, change the dressing. Then press "On" to restart the therapy.</p>
	<p>Recharge the battery either by placing the Invia Liberty pump in the docking station or plug in the charger to the electrical outlet port on the pump.</p> <p>Remaining time of battery is 15 minutes.</p>
	<p>Change the canister, see chapter “Change Invia Liberty canister and Invia Liberty tubing”.</p>

Alarm table

	Fault number	Problem description on the display	Troubleshooting on the display	Pressure
Alarm	311	 <p>Selftest failed</p> <p>311</p> <p>Mute</p>	 <p>Snap the canister out and in again</p> <p>311</p> <p>On</p>	X
	312	 <p>Pump in Standby</p> <p>Leak in system</p> <p>312</p> <p>Mute</p> <p>On</p>	 <p>Check dressing for air leakage and if canister is properly inserted.</p> <p>Consult IFU for further instructions.</p> <p>312</p> <p>On</p>	X
	313	 <p>Filter clogged</p> <p>313</p> <p>Mute</p> <p>On</p>	 <p>Change canister</p> <p>313</p> <p>On</p>	X
	315	 <p>Acceptable internal temperature exceeded</p> <p>315</p> <p>Mute</p>	 <p>Switch pump off and on.</p> <p>If problem persists, contact your Health Care representative.</p> <p>315</p>	X
		 <p>Internal fault</p> <p>316</p> <p>Mute</p>	 <p>Switch pump off and on.</p> <p>If problem persists, contact your Health Care representative.</p> <p>316</p>	X

CAUTION

If fault repeats, note the fault number, switch off the pump and contact the Medela customer service.

CAUTION

Invia Liberty instructions advise 24 hours therapy without interruption. If therapy is discontinued for more than 2 hours using foam or gauze, the dressing should be replaced and therapy restarted by a healthcare professional.

	Remarks/potential cause of fault
	<p>Canister: - Release the canister and reposition. Tubing: - Ensure that the tubing is not twisted, kinked or clamped.</p>
	<p>Occurs when alarm 301 was pending for 5 minutes. For troubleshooting, refer to instructions for Alarm 301. Restart the pump by pressing “On” [].</p>
	<p>To restart the therapy after canister has been changed, press right selection button “On” [].</p>
	<p>Occurs when Warning 406 was pending for 30 minutes. Cool Invia Liberty pump down.</p>
	<p>Restart the pump. If internal fault remains, turn off by pressing [ > 3 seconds] and contact the Medela customer service.</p>

 OPERATION HALTS


 OPERATION CONTINUES

Accessories overview

CAUTION

The Invia Liberty pump is verified within the scope of conformity evaluation and is only to be used with products included in the Invia Liberty NPWT system and distributed by Medela. Medela can only guarantee performance of the system with these products.

*Product not available in all countries

Wound dressings

 Wound dressings to be applied and changed by healthcare professionals only.

The Invia Liberty NPWT system is intended to be used in conjunction with the Invia dressings only. For specific dressing indications, contraindications, warnings and cautions, consult the appropriate Invia dressing instruction for use.

Sterility and requirements for usage

The Invia Liberty disposables products are sterile and single use devices.

Canister tubing		
Canister		
Y-connector		
Drain adapter and drain		
Invia FitPad (tubing)		These are single use devices that should be disposed of after use. If reused performance of the product may deteriorate, cross-contamination may occur.
Wound dressing		

CAUTION

Do not re-sterilise. Do not use if inner package is damaged or opened prior to use.
Do not reuse. If reused performance of the product may deteriorate, cross-contamination may occur.

Cleaning and disinfection

CAUTION

Invia Liberty pump with associated products (docking station, rail holder and charger) should be cleaned / disinfected after every use. Before cleaning the device, unplug the pump from the wall outlet.

	Cleaning	Disinfection	Sterilization	Dishwasher
Pump	✓	✓	✗	✗
Docking station				
Rail holder				
Charger				
Carrying case*	Wipe off with a damp cloth.	Wipe off with a disinfecting agent	Sterilization / cleaning in a dishwasher is NOT allowed.	

* If needed, the carrying case can be washed in washing machine

CAUTION

Invia Liberty pump, docking station, rail holder, charger and carrying case cannot be sterilized. Immersion disinfection, thermal disinfection and ultrasound cleaning are not permitted.

Disinfection (Pump housing, docking station and charger)

CAUTION

Invia Liberty pump can be disinfected with "alcohol".

Disinfection

Invia Liberty can be disinfected with the disinfecting agent group "alcohol".

Do not use other cleaning agents (e.g. Terralin) as they can damage the plastic housing.

Immersion disinfection, thermal disinfection and ultrasound cleaning are not permitted.

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CAUTION

Do not use other cleaning agents (e.g. Terralin) as they can damage the plastic material.

Cleaning procedure for Invia Liberty pump, docking station and charger

1. Wear suitable protection (clothing, gloves, face mask and goggles) according to local guidelines.
2. Apply the disinfectant agent in accordance with instructions from the manufacturer¹⁾. Pay attention particularly to edges, narrow corners and bottom side.
3. Leave the disinfectant on. Follow the recommended residence time for the disinfectant as instructed by the manufacturer¹⁾.
4. Thoroughly clean the surface, all edges, housing niches, corners, bed holder, brackets, port covers and bottom side.

5. Wipe dry or air dry as instructed from the manufacturer¹⁾.
6. If needed, repeat step 2-5 to ensure proper cleaning.
7. Dispose contaminated material in accordance with local environmental guidelines.

**CAUTION**

Do not spray disinfectants directly into openings as this may harm electronic components.

For a detailed cleaning instruction, contact your Medela representative.

¹⁾ Manufacturer of the disinfectant agent.

Disposal

Invia Liberty pump is made from various metals and plastics. Before disposal the rechargeable battery and electronics must be removed according to instructions. Then Invia Liberty pump is no longer operational. Disposal of electronics and plastic components should be handled in accordance with local environmental guidelines.

Invia Liberty disposables and dressings should be handled and disposed in accordance with local environmental guidelines.

**CAUTION**

Must not be disposed of together with household refuse.

Maintenance / Safety-related check

Service work may only be carried out by authorized personnel. A safety-related check is confined to visual inspection of the housing and charger for damage and must be performed prior to each use. If Invia Liberty pump has not been in use, the battery must be charged approximately once every 6 months to ensure optimum function.

Guarantee

Guarantee for 2 years after date of delivery in used in accordance with these instructions. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate usage as well as use by unauthorized persons.

Service life

The service life of the device is five years; the internal batteries life included.

Signs and symbols

This symbol indicates a safety related tip.

This symbol indicates a CAUTION or WARNING associated with the device.

This symbol indicates a class II device.

This symbol indicates to not dispose the device together with household refuse (for EU only).

This symbol indicates the date of manufacture (four digits for the year and two digits for the month).

This symbol indicates the name and the address of the manufacturer.

This symbol indicates the device is sterilized using ethylene oxide.

This symbol indicates a prescription device. CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician. (for US only).

This symbol indicates a type BF applied part.

This symbol indicates manufacturer's catalog number.

This symbol indicates manufacturer's serial number.

This symbol indicates manufacturer's batch code.

This symbol indicates the protection against the ingress of solid foreign objects and against harmful effects due to the ingress of water.

This symbol indicates that the device should not be used after the end of the year and month shown.

This symbol indicates to follow the instruction for use.

This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.

This symbol indicates to not use the device if package is damaged.

This symbol indicates the number of items.

This symbol indicates the direct current socket.

This symbol indicates a single use device. Do not reuse the device.

This symbol indicates MR unsafe.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation and storage.

This symbol indicates the atmospheric pressure limitation for operation and storage.

This symbol indicates to keep the device dry.

This symbol indicates to handle the fragile device with care.

This symbol indicates to keep the device away from sunlight.

This symbol indicates that the device is in conformance with the Medical Device Directive 93/42/EEC.

This symbol indicates that the marked item or its material is part of a recovery or recycling process.

This symbol indicates the correct upright position of the transport package.

Technical specifications

vacuum range
- 40 to -200 mmHg
- 5 to -27 kPa

low flow
5 L/min
flow

without canister
1000 g
2.2 lbs

IP33

H x W x D
150 x 170 x 95 mm


AC ISO 13485
CE (93/42/EEC), IIa


max. noise level
- 42.5 dB(A) 1 l

alarm noise level
78 dB(A)

Operating conditions

Transport/
Storage conditions

Switching adapter AC
Model: MSA-C2500IS12.0-30C-ZZ
IEC: 60601-1
Input: 100-240VAC, max. 0.8A,
50/60Hz
Output: 12.0VDC, 2.5A

Electromagnetic compatibility (EMC)

Invia Liberty is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014 4th Edition according clause 7 and 8.9.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition and IEC 60601-1-2:2014 4th Edition, Table 1)

Electromagnetic emissions		
Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Invia Liberty pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	The Invia Liberty pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

WARNING

The Invia Liberty pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Invia Liberty pump should be observed to verify normal operation in the configuration in which it will be used.

WARNING

Invia Liberty is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information.

WARNING

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Invia Liberty pump.

WARNING

EMC

HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Invia Liberty pump.

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 2)

Electromagnetic immunity			
The Invia Liberty pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Liberty pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line-to-earth	± 1 kV differential mode ± 2 kV line-to-earth not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Invia Liberty pump requires continued operation during power mains interruptions, it is recommended that the Invia Liberty pump be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 4)

Electromagnetic immunity			
The Invia Liberty pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Invia Liberty pump should assure that it is used in such environment.			
Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Invia Liberty pump including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 0.35 \sqrt{\frac{P}{\text{Where } P \text{ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and } d \text{ is the recommended separation distance in meters (m).}}}$ $d = 0.35 \sqrt{\frac{P}{80 \text{ MHz to } 800 \text{ MHz}}}$ $d = 0.7 \sqrt{\frac{P}{800 \text{ MHz to } 2.5 \text{ GHz}}}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Invia Liberty pump is used exceeds the applicable RF compliance level above, the Invia Liberty pump should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating the Invia Liberty pump.</p> <p>^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.</p>			

Electromagnetic compatibility (EMC, IEC 60601-1-2:2007 3rd Edition, Table 6)

Recommended separation distance between portable and mobile RF communications equipment and the Invia Liberty pump				
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter			
	m			
W	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$	
0.01	0.035	0.035	0.07	
0.1	0.11	0.11	0.22	
1	0.35	0.35	0.7	
10	1.11	1.11	2.21	
100	3.5	3.5	7.0	
For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

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Medela AG
Medela UK Ltd.