

Guidelines for Cleaning, Disinfecting
and Sterilizing Medela Healthcare Products

REPROCESSING GUIDE

Contents

Reprocessing Guide (Cleaning/Disinfection/Sterilization)	3
1 General Information.....	3
2 Requirements for Use and Re-Use	3
3 Staff	4
4 Definitions	4
Processing Instructions for Non-Critical Medical Products	
at Health Care Facilities	6
5 List of Non-Critical Medical Products.....	6
6 Processing Instructions.....	9
7 Point of Use.....	10
8 Thorough Cleaning.....	11
9 Intermediate Level Disinfection.....	12
10 Storage	13
Processing Instructions for thermal sensitive Non-Critical Medical Products	
at Health Care Facilities	14
11 List of thermal sensitive Non-Critical Medical Products.....	14
12 Processing Instructions.....	15
13 Point of Use.....	16
14 Thorough Cleaning and Thermal Disinfection	16
15 Storage	18
Processing Instructions for re-sterilizable Semi-Critical Medical Products	
at Health Care Facilities	19
16 List of re-sterilizable Non-Critical Medical Products.....	19
17 Processing Instructions.....	21
18 Point of Use.....	22
19 Thorough Cleaning and Thermal Disinfection	23
20 Terminal Sterilization.....	25
21 Storage	26
Disassembly and Reassembly Chart	27
References	27
Phone number	28

Reprocessing Guide (Cleaning / Disinfection / Sterilization)


1 General Information

This Guide is intended to give general instructions on how medical devices supplied by Medela AG may be processed or reprocessed to prepare them for Re-Use and gives advice on cleaning, disinfecting or sterilizing products supplied by Medela AG.

Equipment, operators, cleaning detergents / disinfectant agents and procedures have a contribution to the efficacy of the processing and reprocessing of medical products.

These guidelines are for healthcare professionals, such as nurses and infection control practitioners, who work in a hospital environment.

2 Requirements for Single Use and Re-Use

Type	Use for	Requirements
Re-Use <i>e.g.</i> <i>tubing</i>	↻ more than one patient multiple uses possible	<ul style="list-style-type: none">- Requires reprocessing after each use- Must be checked regularly for damage and appearance and replaced if necessary
Single Use  <i>e.g.</i> <i>filters</i>	↻ one patient only once only	<ul style="list-style-type: none">- Ready for use- Use immediately after opening the packaging- Must be disposed after using- Not to be used by more than one person in order to avoid possible health risks- Reuse could cause cross contamination- Reprocessing could cause loss of mechanical, chemical and / or biological characteristics

Note:

Medela AG does not define the maximum number of uses appropriate for reprocessing and reusing of medical devices. The useful life of these devices depends on many factors including the method and duration of each use, and the handling between uses. Careful inspection and functional test of the device before use is the best method of determining the end of serviceable life for the medical device. Refer to Chapter 19 of the suction pumps' Instructions for Use for additional re-placement guidelines.

Not all products are available in all markets.

3 Staff

Staff should use suitable personal protective clothing and equipment at all times. In particular take note of the instructions provided by the cleaning detergent / disinfecting agent manufacturer for correct handling and use of the product.

4 Definitions

4.1 Material abbreviations

ABS	Poly(Acrylonitrile Butadiene Styrene)
ASA	Poly(Acrylic Styrene Acrylonitrile)
PA	Polyamide
PC	Polycarbonate
POM	Polyoxymethylene (Acetal)
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
TPE	Thermoplastic Elastomer

4.2 Point of Use

Reprocessing already begins at the point of use, which includes initial cleaning and measures to prevent drying of soil and contaminants in and on the device and should be carried out directly after application (within a maximum of 2 hours after use).

4.3 Initial Cleaning and Rinsing with sterile cold Water

The first step of decontamination is initial cleaning with sterile cold water (< 40 °C, < 104 °F). The purpose of the initial cleaning is the removal of all visible soils from an item and the reduction of the number of particulates, microorganisms and potential pyrogens. Many soils, including proteinaceous blood components such as albumin and hemoglobin, are water-soluble and can be easily washed away with water alone.

4.4 Thorough Cleaning

Thorough cleaning consists of the removal, usually with detergent and water, of adherent soil (e.g., blood, protein substances and other debris) from the surfaces, crevices, serrations, joints, and lumens of instruments, devices and equipment by a manual or mechanical process that prepares the items for safe handling and/or further decontamination (Disinfection / Terminal Sterilization).

4.5 Disinfection

Disinfection is a process that removes, inactivates or destroys blood-borne pathogens and other microorganisms by physical or chemical means on a surface or item to the point where they are no longer capable of transmitting infectious particles and the surface or item is rendered safe for handling, for use or for disposal. Disinfection processes do not ensure the same margin of safety associated with sterilization processes.

4.6 Thorough Disinfection or Terminal Sterilization

The complete elimination or destruction of all forms of microbial life accomplished in health care facilities by either physical or chemical processes.

Processing Instructions for Non-Critical Medical Products at Health Care Facilities

5 List of Non-Critical Medical Products

Recommended cleaning methods have been validated by Medela AG to current international guidelines for cleaning.

Electric and electronic equipment and accessories

Picture	Item No	Description	Material
	014.0XXX	Clario / Clario Toni Motor unit AC-version Motor unit AC/DC-version	Housing: Plastics (ABS)
	074.0006	Clario / Clario Toni vacuum gauge	Metals (stainless) Plastics (various)
	074.0007	Clario / Clario Toni carrying bag	Plastics (PE coated with PET)
	026.XXXX	Vario 8 AC + AC/DC Vario 8 c/i AC/DC Vario 18 AC + AC/DC Vario 18 c/i AC/DC	Housing Plastics (ABS)
	077.0105	Vario jar adapter	Plastics (ABS)
	077.0511	Trolley small	Metals (Zinc-coated steel)
	077.0821	Vario rail holder	Metals (AlMgSi)
	077.0823	Vario universal holder	Metals (Aluminum)

Picture	Item No	Description	Material
	077.0831 077.0832	Vario carrying bag Vario car connection cable	Plastics (Polyester with PVC coating)
	087.0000	Invia Liberty pump	Housing: Plastics (ABS) Holder: Plastics (PAA) Button: Plastics (POM)
	079.0037	Docking station	Housing: Plastics (PC) Knob: Plastics (PA66)
	077.0148	Mains adapter winternational	n/a
	071.0000 071.0001	Basic suction pump, rack and portable versions	Housing Plastics (ABS)
	071.0002 071.0003	Dominant Flex suction pump, rack and portable versions	Housing Plastic (ABS)
	071.0034 071.0035	Trolley to Dominant Flex and Basic suction pumps	PS-coated metal
	077.0104 077.0521 077.0152	Clampholder	Plastics (PA) Metals (Aluminum)
	077.0723	Foot On/Off switch	Various plastics

Picture	Item No	Description	Material
	077.0731	Foot vacuum regulator	Various plastics
	077.0751 077.0752	Tubing holders	Metals (stainless)
	017.0015	Apgar Timer	Housing: Plastics (PC/ABS)
	077.0523	Wall bracket, single	Plastics (ASA)
	079.0000 079.0002 079.0003	Thopaz drainage pump	Housing: Plastics (ABS)
	079.1000 079.1002 079.1003	Thopaz+ digital cardio-thoracic drainage system	Housing: Plastics (ABS)
	079.0031	Carrying strap Thopaz	Textiles, Plastics (ABS)
	079.1006	Carrying strap Thopaz+	Textiles, Plastics (ABS)
	079.0036 079.0040	Holder with standard rail	Metals (stainless, aluminium)
	079.0038	Adapter docking station for Thopaz canister 2l	Plastics (PP)
	077.1456	Vacuum gauge	Metals (stainless) Various plastics

6 Processing Instructions

The sequence of steps required to prepare medical devices for reprocessing are summarized in the chart below.

More detailed instructions for each step are given on the following pages.

Point of Use

- Disassembly
- Manual Initial Cleaning



Thorough Cleaning

- Manual Thorough Cleaning
- Drying




Disinfection

- Manual Disinfection
- Drying



Storage

- Inspection
- Reassembly



7 Point of Use

7.1 Utensils and Equipment

- Personal Protective Equipment (PPE: Disposable gloves, Proper protective gear).
- Lint free nonwoven wipes.
- Sterile cold water (< 40 °C, < 104 °F).

7.2 Procedure Disassembly

- Disinfect your hands.
- Put on disposable gloves and proper protective gear.
- Disassemble if possible.



Housing suction pumps: remove Safety Set as shown.
Suction pumps in general: Remove all accessories if attached.

7.3 Procedure Manual Initial Cleaning

- Directly after application (within a maximum of 2 hours postoperatively) wipe the surfaces to remove all gross debris by using a lint free nonwoven wipe wetted with sterile cold water (< 40 °C, < 104 °F).
- Ensure that all surfaces are thoroughly wetted.
- Wipe again using another wipe and keep moist for a minimum of 5 minutes then repeat it with another wipe.
- Dispose of protective gloves.
- Disinfect your hands.

7.4 Cautions

- Before cleaning electrical devices, disconnect the power plug from the power socket.
- Clean all surfaces immediately after use (maximum 2 hours postoperatively) to avoid residues from drying and to prevent growth of microorganisms.
- Never use steel brushes or steel wool for cleaning.
- Water temperatures higher than 45 °C (113 °F) can cause protein coagulation which may affect the efficacy of the process.
- Never place the device in water or other liquids.

8 Thorough Cleaning

8.1 Utensils and Equipment

- Personal Protective Equipment (PPE: Disposable gloves, Proper protective gear).
- Lint free nonwoven wipes wetted with an enzymatic detergent.
- Sterile cold water (< 40 °C, < 104 °F).

8.2 Procedure Manual Thorough Cleaning

- Disinfect your hands.
- Put on disposable gloves and proper protective equipment.
- Use the wipes wetted with an enzymatic detergent according to the manufacturer's instructions for use.
- Ensure that all surfaces are thoroughly wetted.
- Check surfaces for visible dirt and repeat these steps if necessary.
- At the end use a wipe wetted with sterile cold water (< 40 °C, < 104 °F) to wipe surfaces. Repeat these steps if necessary.

8.3 Procedure Drying

- Allow surfaces to dry from the enzymatic detergent according to the manufacturer's instructions for use.
- Dispose of protective gloves.
- Disinfect your hands.

8.4 Cautions

- Before cleaning electrical devices, disconnect the power plug from the power socket.
- Never use steel brushes or steel wool for cleaning.
- Water temperatures higher than 45 °C (113 °F) can cause protein coagulation which may affect the efficacy of the process.
- Never place the device in water or other liquids.

9 Intermediate Level Disinfection

9.1 Utensils and Equipment

- Personal Protective Equipment (PPE: Disposable gloves, Proper protective gear).
- Disinfecting wipes.

9.2 Procedure Manual Disinfection

- Disinfect your hands.
- Put on disposable gloves and proper protective equipment.
- Use disinfecting wipes according to the manufacturer's instructions for use.
- After 5 minutes wipe with another disinfecting wipe.

9.3 Procedure Drying

- Allow the surface to dry for a minimum of 5 minutes.
- Dispose of protective gloves.
- Disinfect your hands.

9.4 Caution

- Never place the device in water or other liquids.

9.5 Recommended agent:

- CaviWipes®
Metrex® Research
Address: 1717 West Collins Avenue, Orange, CA 92867, U.S.A.
Homepage: <http://www.metrex.com>
Phone: (800) 841 1428
Email: metrexcustcare@sybrondental.com
- Mikrozyd® AF Wipes
Schülke&Mayr GmbH
Address: Robert-Koch Str. 2, 22851 Norderstedt, GERMANY
Homepage: <http://www.schuelke.com>
Phone: +49 (0) 40 521 00 0
Email: info@schuelke.com

10 Storage

10.1 Procedure Reassembly and Storage

- Where appropriate the cleaned, disinfected and checked medical products should be reassembled again and prepared for Re-Use.



Housing suction pumps: Reattach Safety Set as shown
Suction pumps in general: Configure pump with necessary accessories

- Sufficient protection of the medical products must be provided in order to keep medical products dry and dust free.

10.2 Cautions






- Store medical products dry and dust free.

Processing Instructions for thermal sensitive Non-Critical Medical Products at Health Care Facilities

11 List of thermal sensitive Non-Critical Medical Products

Recommended cleaning methods have been validated by Medela AG to current international guidelines for cleaning.

Thermo sensitive medical products

Picture	Item No	Description	Material
	074.0001	Clario/Clario Toni canister Set consists of canister, lid with plugs for closing and floater (074.0010)	Jar: Plastics (PC) Lid: Plastics (TPE-E) Plugs: Plastics (PE) Floater: Plastics (TPE-V)
	074.0002 014.0030	Clario/Clario Toni SafetyChamber is a set consisting of lid, housing, 2 valves (074.0008)	Housing: Plastics (PBT+ASA) Lid: Plastics (TPE-E) Valves: Elastomers (MQ)
	014.0005	Clario/Clario Toni vacuum regulator Knob	Plastics (TPE-V)
	077.0082 077.0085	PC suction jar, 1.5l PC suction jar, 2.5l	Plastics (PC)
	077.0531 077.0532	Quiver, 280 ml with holder Quiver, 480 ml with holder	Plastics (PP)

12 Processing Instructions

The sequence of steps required to prepare medical devices for reprocessing are summarized in the chart below.

More detailed instructions for each step are given on the following pages.

Point of Use

- Disassembly
- Manual Initial Cleaning



Thorough Cleaning

- Automated Cleaning with Washer Disinfector



Thermal Disinfection

- Automated Disinfection with Washer Disinfector
- Inspection



Storage

- Reassembly
- Packaging
- Storage



13 Point of Use

13.1 Utensils and Equipment

- Personal Protective Equipment (PPE: Disposable gloves, Proper protective gear)
- Sterile cold water (< 40 °C, < 104 °F)

13.2 Procedure Disassembly

- Disinfect your hands.
- Put on disposable gloves and proper protective equipment.
- PC jars:
 - When used with the disposable liners from Medela: remove liners according to suction pump's Instructions for Use, chapter 10
 - Remove jar from clamp holder
- Quiver:
 - Remove from holder

13.3 Procedure Initial Cleaning with Rinsing Water

- Remove gross debris.
- Rinse all parts with sterile cold water (< 40 °C, < 104 °F) for 3 minutes.
- Dispose of protective gloves.
- Disinfect your hands.

13.4 Cautions

- Clean all parts immediately after use (maximum 2 hours postoperatively) to avoid drying of residues and to prevent growth of microorganisms.
- Never use steel brushes or steel wool for cleaning.
- Water temperatures higher than 45 °C (113 °F) can cause protein coagulation which may affect the efficacy of the process.

14 Thorough Cleaning and Thermal Disinfection

14.1 Utensils and Equipment

- Disposable gloves

14.2 Procedure Automated Cleaning and Disinfection with Washer Disinfector

- Disinfect your hands.
- Put on disposable gloves.
- Load items into the washer disinfector.
- Choose the washer disinfector washing cycle.
- Operate the washer disinfector.

- On completion, unload the washer disinfectant.
- Visually check for remaining soil and dryness. If soil remains, repeat the procedure Cleaning and Disinfection with Washer Disinfectant.
- Remaining wetness may be removed with medical grade compressed air or by heating in an oven below 110 °C (230 °F).

14.3 Procedure Inspection

- Check the parts visually for damage and throw away at first signs of damage or weakness.
- Dispose of protective gloves.
- Disinfect your hands.

14.4 Cautions

- Water temperatures higher than 45 °C (113 °F) during the cleaning cycle can cause protein coagulation which may affect the efficacy of the process.
- Avoid contact between items (movement during washing could cause damage).
- Arrange items so that channels and openings are always oriented downwards
- Please follow Instructions for Use of manufacturer:
 - for washing disinfectant and
 - for detergents and disinfectants
- If available, attach tubing to nozzles in washer disinfectant.

14.5 Recommended agents for thorough cleaning and thermal disinfection:

- neodisher® FA as the Liquid Detergent and neodisher® Z as the Liquid Neutralizing Agent Miele, Inc.

Address: 9 Independence Way, Princeton, NJ 08540, U.S.A.

Homepage: <http://www.miele.us>

Phone: (800) 991 9380

Email: proinfo@mieleusa.com

- neodisher® FA as the Liquid Detergent and neodisher® Z as the Liquid Neutralizing Agent Chemische Fabrik Dr. Weigert GmbH & Co. KG Mühlenhagen 85, D-20539 Hamburg, GERMANY
Homepage: <http://www.drweigert.com>

15 Storage

15.1 Procedure Reassembly and Storage

- Where appropriate, the cleaned, disinfected and checked medical products should be prepared for Re-Use.
- Sufficient protection of the medical products must be provided in order to keep medical products dry and dust free.

15.2 Cautions

- Store medical products dry and dust free.









Processing Instructions for re-sterilizable Non-Critical Medical Products at Health Care Facilities

16 List of re-sterilizable Non-Critical Medical Products

Recommended cleaning methods have been validated by Medela AG to current international guidelines for cleaning.

Sterilizable Medical Products

Picture	Item No	Description	Material
	077.0420 / 077.0430 077.0440 / 077.0450	Lid	Plastics (PA, TPE, POM)
	077.0110 / 077.0120 077.0125 / 077.0130 077.0150 / 077.0155	Suction jar PSU	Plastics (PSU)
	077.0900 / 077.0901 077.0902 / 077.0905 077.0970 / 077.0095 077.0542 / 077.0912 077.0913 / 077.0921 077.0922 / 077.0931 077.0961 / 077.0185	Silicone tubing	Elastomers (Silicone) Plastics (POM) Elastomers (silicone o-ring)
	077.1018 077.1019	Coupling pieces	Plastics (POM) Elastomers (silicone o-ring)
	077.1022	Drainage valve	Plastics (PSU) Elastomers (silicone o-rings, bellows)
	077.0563	Change- over valve	Plastics (PSU, PA) Elastomers (silicone tubing)

	<p>077.0651 / 077.0701 077.0702 / 077.0703 077.0704 / 077.0705 077.0706 / 077.0707 077.0715 / 077.0716 077.0711</p>	<p>Assembled products and sets</p>	<p>Jar: PSU Lid: PA, TPE, POM</p>
	<p>077.0029 / 077.0030 077.0031 / 077.0038 077.0175 / 077.0176 077.0177</p>	<p>Malmström Cups</p>	<p>Metals (stainless) Elastomers (silicone)</p>
	<p>077.1041 / 077.1042 077.1043 / 077.1044 077.0173 / 077.0172 077.0171</p>	<p>Bird Cups</p>	<p>Metals (stainless) Plastics (PP, POM) Tubing: Elastomeres (Silicone)</p>
	<p>077.0063 077.0174</p>	<p>Bird Posterior Cup</p>	<p>Metals (stainless) Plastics (PP, POM) Tubing: Elastomeres (Silicone)</p>
	<p>077.0078</p>	<p>Reusable Silc Cup with trump- pet valve</p>	<p>Metals (stainless) Elastomers (silicone)</p>
	<p>077.0761 077.0762</p>	<p>Reusable Silc Cups</p>	<p>Metals (stainless) Elastomers (silicone)</p>
	<p>077.0080</p>	<p>CaesAid Cups for Caesarean sections</p>	<p>Elastomers (silicone)</p>
	<p>077.0081</p>	<p>Manual Vacuum Extractor</p>	<p>Metals Plastics (PP) Elastomers (silicone)</p>

17 Processing Instructions

The sequence of steps required to prepare medical devices for reprocessing are summarized in the chart below.

More detailed instructions for each step are given on the following pages.

Point of Use

- Disassembly
- Manual Initial Cleaning



Thorough Cleaning

- Automated Cleaning with Washer Disinfector




Thermal Disinfection / Sterilization

- Automated Disinfection with Washer Disinfector
- Inspection, Packaging
- Terminal Sterilisation by Steam



Storage

- Storage



18 Point of Use

18.1 Utensils and Equipment

- Personal Protective Equipment (PPE: Disposable gloves, Proper protective gear)
- Sterile cold water (< 40 °C, < 104 °F)

18.2 Procedure Disassembly

- Disinfect your hands.
- Put on disposable gloves and proper protective equipment.
- Disassemble the product (according to the instructions sheets delivered with the product) into individual parts.
- Change-over valve:
 - Disassemble according to instruction sheet supplied with product

18.3 Procedure Initial Cleaning with Rinsing Water

- Rinse all parts with sterile cold water (< 40 °C, < 104 °F) for 3 minutes
- Dispose of protective gloves.
- Disinfect your hands.

18.4 Cautions

- Clean all parts immediately after use (maximum 2 hours postoperatively) to avoid drying of residues and to prevent growth of microorganisms.
- Never use steel brushes or steel wool for cleaning.
- Water temperatures higher than 45 °C (113 °F) can cause protein coagulation which may affect the efficacy of the process.

19 Thorough Cleaning and Thermal Disinfection

19.1 Utensils and Equipment

- Disposable gloves

19.2 Procedure Automated Cleaning and Disinfection with Washer Disinfector

- Disinfect your hands.
- Put on disposable gloves.
- Load items into the washer disinfector.
- Choose the washer disinfector washing cycle.
- Operate the washer disinfector.
- On completion unload the washer disinfector.
- Visually check for remaining soil and dryness. If soil remains repeat the procedure Cleaning and Disinfection with Washer Disinfector.
- Remaining wetness may be removed with medical grade compressed air or by heating in an oven below 110 °C (230 °F).

19.3 Procedure Inspection

- Check the parts visually for damage and throw away at first signs of damage or weakness.
- Dispose of protective gloves.
- Disinfect your hands.

19.4 Cautions

- Water temperatures higher than 45 °C (113 °F) during the cleaning cycle can cause protein coagulation which may affect the efficacy of the process.
- Avoid contact between items (movement during washing could cause damage)
- Arrange items so that channels and openings are always oriented downwards.
- Please follow Instructions for Use of manufacturer:
 - for washing disinfector and
 - for detergents and disinfectants
- If available, attach tubing to nozzles in washer disinfector.

19.5 Recommended agents for thorough cleaning and thermal disinfection:

- neodisher® FA as the Liquid Detergent and
neodisher® Z as the Liquid Neutralizing Agent

Miele, Inc.

Address: 9 Independence Way, Princeton, NJ 08540, U.S.A.

Homepage: <http://www.miele.us>

Phone: (800) 991 9380

Email: proinfo@mieleusa.com

- neodisher® FA as the Liquid Detergent and
neodisher® Z as the Liquid Neutralizing Agent

Chemische Fabrik Dr. Weigert GmbH & Co. KG

Address: Mühlenhagen 85, D-20539 Hamburg, GERMANY

Homepage: <http://www.drweigert.com>

Phone: +49 (0) 407 89 60 0

Email: info@drweigert.de

20 Terminal Sterilisation

20.1 Utensils and Equipment

- Steam Steriliser
- Sterilisation cycle Information

Method	Moist Heat Sterilisation according to ANSI/AAMI ST79
Cycle	Gravity-displacement Steam Steriliser
Temperature	135 °C (275 °F)
Exposure Time	10 minutes
Pressure	0.1 to 2.2 bar (2 to 22 psia)
Drying Time	20 minutes (minimum, in chamber)
Cooling Time	20 minutes (minimum, at room temperature)

20.2 Procedure Terminal Sterilisation

- Visually check for remaining soil and dryness. If soil remains repeat the procedure Cleaning and Disinfection with Washer Disinfector.
- Handle items that need terminal sterilizing according to internal hospital guidelines.
- Load items into the Steam Steriliser.
- Choose the sterilisation cycle.
- Operate the Steam Steriliser
- On completion unload the Steam Steriliser.

20.3 Cautions

- Do not stack items in the steriliser. Place items correctly and loosely into the basket, shelf, or cart of the steriliser so as not to impede the penetration of the sterilant.
- Please follow Instructions for Use of manufacturer for Steam Steriliser.

21 Storage

21.1 Procedure Storage

- Where appropriate the cleaned, disinfected and checked medical products should be assembled for re-use according to the instruction sheets delivered with assembled products.
- Sufficient protection of the medical products must be provided in order to keep medical products dry and dust free.
- Ensure that the sterile storage area is a well-ventilated area that provides protection against dust, moisture, insects, and temperature and humidity extremes.
- Store sterile items so the packaging is not compromised (e.g. punctured, bent). After disinfection, store the components in clean foil until required for use.
- Label sterilized items with a load number that indicates the steriliser used, the cycle or load number, the date of sterilisation, and, if applicable, the expiration date according to in-house hospital policy.
- The shelf life of a packaged sterile item depends on the quality of the wrapper, the storage conditions and the conditions during transport, the amount of handling, and other events (moisture) that compromise the integrity of the package.
- Evaluate packages before use for loss of integrity (e.g., torn, wet and punctured). The pack can be used unless the integrity of the packaging is compromised, repack and reprocess before use.

21.2 Cautions

- Store medical products dry and dust free.

Disassembly and Reassembly Chart

See instructions from the suction pump's instructions for use and individual instruction sheets supplied with accessories for disassembly.

References

- U.S. Food and Drug Administration (FDA)
Draft Guidance for Industry and FDA Staff: May 2, 2011
Processing/Reprocessing Medical Devices in Health Care Settings:
Validation Methods and Labeling
- AAMI TIR12:2004, Arlington, VA: AAMI, 2004. Technical Information Report
Association for the Advancement of Medical Instrumentation
Designing, testing, and labeling reusable medical devices for reprocessing in
health care facilities: A guide for medical device manufacturers
- AAMI TIR30:2003, Arlington, VA: AAMI, 2003. Technical Information Report
Association for the Advancement of Medical Instrumentation
A compendium of processes, materials, test methods, and acceptance
criteria for cleaning reusable medical devices
- ANSI/AAMI ST79: 2010 & A1: 2010 & A2:2011 & A3:2012
Association for the Advancement of Medical Instrumentation
Comprehensive Guide to Steam Sterilisation and Sterility Assurance in Health
Care Facilities
- Centers for Disease Control and Prevention (CDC)
Guideline for Disinfection and Sterilisation in Healthcare Facilities, 2008
- EN ISO 17664:2004:
Sterilisation von Medizinprodukten
Vom Hersteller bereitzustellende Informationen für die Aufbereitung von
resterilisierbaren Medizinprodukten

Phone numbers

For additional information on the device, including questions on infection control procedures, please call:
USA (+1) 877 735 1626, International (+41) 41 562 51 51



Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
www.medela.com

International Sales

Medela AG
Lättichstrasse 4b
6341 Baar, Switzerland
Phone +41 41 562 51 51
Fax +41 41 562 51 00
ism@medela.ch
www.medela.com

Local contact: